


## No Complement Fixation Test to detect antibodies against *Chlamydia* any longer

Unfortunately, the manufacturer of CFT reagents has now stopped producing essential reagents. Therefore, from March 15, 2025, we will no longer be able to carry out the complement fixation reaction (CFT) to detect antibodies against *Chlamydia*.

**For the detection of infections in pigs, ruminants and also pets (dogs and cats), we offer direct detection of *Chlamydia species* using real-time PCR.**

Depending on the clinical symptoms (reproductive disorders or conjunctivitis), specimens should be taken from the affected organ system (genital swab, conjunctival swab, placenta, fetal liver, lung, spleen, genital organs).

If you have any further questions about *Chlamydia* diagnostics, please call us on: +49-511-220029-0



IVD Gesellschaft für Innovative Veterinärdiagnostik mbH  
Albert-Einstein-Straße 5  
D- 30926 Seelze  
Telefon: +49-511-220029-0  
service@ivd-gmbh.de  
www. ivd-gmbh.de

**Status 3/2025**